




The Saviour Tactical is a roll up stretcher that only weighs 5 kg. It is easy to carry, easy to deploy, and easy to re-stow. Every Saviour Tactical stretcher has the following qualities

- Works well in a wet or sandy environment
- Neutrally buoyant - so floats if dropped overboard.
- Will float with a patient in a lifejacket
- Excellent for confined space rescue
- Features cross chest straps and figure of 8 foot harness to allow vertical lifting without fear of patient sliding out of either end of the stretcher
- Has a lifting eye at the top for vertical haul rated at 200kg
- Head and foot end drag harnesses for single man drag rescues
- 6 soft lined hand hold for team carries
- Can be carried by 2 or 4 people, +/- the addition of poles through the handles
- Head hugger system to protect and immobilise head, which can be used with or without a collar.
- Extra central layer of skin material to increase spinal support
- Acts as an effective pelvic splint
- Xray translucent and CT friendly
- Discreet colour coded straps for the tactical environment (although these can be made more obvious for non-tactical environments)
- All straps are replaceable independently
- Easily infection control compliant meaning that it is not limited to single use.

SAVIOUR TACTICAL *stretchers*



Saviour Tactical Stretcher: IDENTIFICATION OF PARTS LIST



Non-replaceable items

- A. MAIN SKIN
- B. INNER SKIN

Replaceable straps:

1. Chest strap – closure side with handle
- 1a. Chest strap – buckle side with handle
2. Waist strap – closure side
- 2a. Waist strap – buckle side
3. Hip strap – closure side with handle
- 3a. Hip strap – buckle side with handle
4. Foot loop – closure side
- 4a. Foot loop – buckle side
5. Leg strap – closure side with handle
- 5a. Leg strap – buckle side with handle
6. Left shoulder strap – closure side
- 6a. Left shoulder strap – buckle side
7. Right shoulder strap – closure side
- 7a. Right shoulder strap – buckle side

Other replaceable items:


8. Head Hugger
9. Chin strap
10. Top drag harness
11. Bottom drag harness
12. Closure strap
13. Inner skin retainer
14. Stretcher carry bag

Construction materials:

- All closure, carry straps handles and stitching are POLYPROPYLENE
- All buckles and keepers are 316 STAINLESS STEEL
- Main body and inner skin are LDPE
- Head Hugger and Carry Bag are PVC
- Head Haul Eyelet is 316 STAINLESS STEEL and NYLON 66
- Eyelets for Foot Loop are NICKLE PLATED BRASS
- Head Haul Strap is POLYPROPYLENE SUPER TAPE
- Inner Skin retainer and screw is NYLON 66

SAVIOUR TACTICAL RANGE

Key Features




Saviour Tactical Stretcher: USE OF THE DRAG HARNESS

There are two drag harnesses, one at the head end and one at the foot end. On removing the protective sheath, a tape loop is released that is permanently attached through metal eyelets on the stretcher.

Placing the loop over the head and around the shoulders will allow the rescuer to drag or pull the casualty in the required direction while having their hands free to carry out other duties.


The same dragging or pulling technique is used when moving the casualty head or feet first.

FORWARD DRAG




direction of travel →

REVERSE PULL




Saviour Tactical Stretcher: CARRYING THE CASUALTY

2 man carry




direction of travel →

4 man carry




direction of travel →

2 (or 4) man pole carry for long distance carrying



2 (or 4) man pole carry for long distance carrying

4 man side to side pole carry for heavy casualties




4 man side to side pole carry for heavy casualties

Saviour Tactical Stretcher: VERTICAL HAUL

The Saviour Tactical stretcher is rated for vertical lifting from the top eyelet (not the drag harness). As this only effects a single point of contact with the casualty, a second safety line secured to the casualty is always recommended.

In order to make the vertical haul easier, it is also recommended that both the upper and lower drag harnesses be removed from the stretcher before the lift is made.



Saviour Tactical Stretcher: IN WATER

Being neutrally buoyant, the Saviour Tactical stretcher floats without assistance.

A casualty can be placed in the Saviour Tactical stretcher with a lifejacket fitted and then swum/towed to the point of rescue using the drag harness. Casualties can also be put into the Saviour Tactical stretcher in the water, (as long as the sea is calm enough), and a lifejacket is fitted first. This would be ideal for spinal casualty recovery in calm water.

